DILUTION INTRUCTIONS & HOW TO USE: High Level = 1:10 General use = 1:50 Fogging = 1:10 - 1:50 Wear suitable protective clothing and apparatus as appropriate, refer to MSDS, Average contact time 2 mins, for details of efficacy data please contact info@steri-7.com

HANDLING: Avoid contact with eyes. Do not ingest.

CAUTION: In case of contact with eyes, rinse immediately with plenty of water and seek medical attention. If swallowed, seek medical attention and show this label.

STORAGE: Keep in original containers. Store in a dry place away from direct sunlight between 4°C and 50°

ACTIVE INGREDIENTS: Benzalkonium Chloride < 1%, Didecyldimethylammonium Chloride <5. PHMB < 1%

## APPROVED TO: ISO 13485, ASTME1052(2000)



For more information please see technical sheet or visit www.steri-7.com or email info@steri-7.com

PCS 95528 CON5





Steri-7 Worldwide Ltd, B1 Aerodrome Studios, 2-8 Airfield Rd. Christchurch BH23 3TH, England. www.steri-7.com

USE DISINFECTANTS SAFELY -Always read the label and product information before use, Never Mix STERI-7 with other products or chemicals.

- ▶ REACTIVE BARRIER TECHNOLOGY - PROTECTS RETWEEN CLEANS
- ► EXCELLENT CLEANING ABILITY
- PROVEN WIDE ANTIMICROBIAL SPECTRUM - BACTERIACIDAL. SPORICIDAL, YEASTICIDAL AND ACTIVE AGAINST ENVELOPED VIRUSES
- ▶ RAPID IN ACTION KILLS UP TO 99.9999% IN SECONDS
- PROVEN ACTIVE IN THE PRESENCE OF ORGANIC MATTER
- DOES NOT DAMAGE MEDICAL DEVICES
- CLEANS AND DISINFECTS AT THE SAME TIME
- ► FULLY EUROPEAN COMPLIANT



GHS05

GHS09

EUH208: Contains polyhexamethylene biguanide hydrochloride. May produce an allergic reaction

H315: Causes skin irritation

H318: Causes serious eye damage. H411: Toxic to aquatic life with long lasting effects.

P280: Wear protective gloves/protective clothing/eye protection/ face protection

P302+P352: IF ON SKIN: Wash with plenty of water. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing P310: Immediately call a POISON CENTER or doctor

P321: Specific treatment (see supplemental first aid instruction on



